

K&A Laborgeraete GbR

Software manual

Synthesizer

H-8


V. 02.11

- I. General
- II. Commands
- III. Settings
- IV. Subroutine programming


I. General

The password is necessary to get into the programming level of the instrument.

- Click 'About' Menu in main Menu -The 'About' Menu will open
- Double Click on the text field 'H - 8' on the 'About' mask (see the drawing)
- 'Password entry' mask will be opened



,About' mask



- Enter the password **2468**

- Click OK button


II. Commands

When password has been entered, new commands in the main menu are accessible:

- 'File – Open - Subroutine' in 'File' submenu
- 'Settings' in main menu




III. Settings



Menu 'Settings' contains 3 submenus which enable entering different settings and parameters which are necessary for instrument control.

- Method
- Directory
- Com

III-1 Settings - Method



This menu enables creating the method used for the synthesis. One method consists of three components: **CYCLE**, **BRANCH TABLE** and **PARAMETERS**.

The single coupling cycle consists of the subroutines which will be executed in the order programmed in the **CYCLE** table.

The Cycle consists of maximum 15 subroutines executed step by step for every synthesis coupling step.

All defined and accessible subroutines are listed in the window on the right side of the **CYCLE**. A subroutine can be chosen by clicking on the name of the desired subroutine.


The synthesis execution starts always from step 1 of the **CYCLE** and continues until the last step of the Cycle. The Cycle will be repeated until the end of the oligo synthesis.

Each subroutine can be repeated 1 to 5 times, according to the entry in the Qty field.


In the example above the subroutine's order is: **1udet 1ucpl 1ucap 1uoxi**.

The **BRANCH TABLE** enables creating the so-called 'Branch'. 'Branch' is the subroutine which can be executed dedicated to the amidite position. The implementing of the 'Branch' will be discussed in chapter **IV-1-4. Subroutine programming-examples** of this manual.

5 different branches can be defined for every amidite position and every wobble combination.




The choice of the amidite position, for which the branch will be programmed, can be made in the choice line on the left top position of the menu.



The choice of the subroutine can be made by clicking on the name of the desired subroutine. All defined and accessible subroutines are listed in the window on the left side of the branch table.

The last component of the Method are **PARAMETERS**




Following parameters are implemented and can be entered:


- **1-Col. Mode:** defines the fluid delivery to the synthesis columns. When YES,

the columns will be filled in serial mode (closer explanation in the chapter **IV-1-1. Subroutine programming-examples**). Status can be changed by clicking at the field.


- **Switch speed AMD:** defines the switching speed of the valves while mixed amidite and other reagent's flow during the synthesis (closer explanation in the chapter **IV-1-2. Subroutine programming-examples**). Lower the number-higher the switch speed of the valves. When 0, all involved valves will be opened and no switching will be executed.
- **Switch speed Reag:** defines the switching speed of the valves while two mixed reagents flow during the synthesis (closer explanation in the chapter **IV-1-2. Subroutine programming-examples**). Lower the number-higher the switch speed of the valves. When 0, all involved valves will be opened and no switching will be executed.
- **Collect DMT:** defines the collecting and evaluation of the trityl monitor during the synthesis. Can be changed to NO by click on the field. Status NO. Recommended only when trityl monitor is not working properly.
- **Add time amd flow:** defines the correction of the flow time when amidite must be delivered to more than two columns at the same time (closer explanation in the chapter **IV-1-3. Subroutine programming-examples**)
- **Add time reag flow:** defines the correction of the flow time when reagent must be delivered to more than two columns at the same time (closer explanation in the chapter **IV-1-3. Subroutine programming-examples**)
- **Add time tet/amd flow:** defines the correction of the flow time when the mixture of amd and other reagent must be delivered to more than two columns at the same time (closer explanation in the chapter **IV-1-3. Subroutine programming-examples**)
- **Add time mix reag flow flow:** defines the correction of the flow time when the mixture of two reagents must be delivered to more than two columns at the same time (closer explanation in the chapter **IV-1-3. Subroutine programming-examples**)
- **AMD dissolve time:** Not applicable in this software version.

When CYCLE, BRANCH TABLE and PARAMETERS are entered, the METHOD can be save now. For that, the name of the Method must be entered in the field: Save as and the key OK must be clicked.






The click on the Methods list enables the selecting of the existing method:




- The click on the button PRINT causes printing of the Cycle, branch table and parameters for selected method.
- The click on the button DELETE causes deleting of the selected method.
- The click on the button EXIT causes exit from the 'METHOD' mask without saving.

III-2 Settings Directory




This menu enables the choice of the folder in which all synthesis data will be managed. The choice can be made by double click on the proper folder and confirming by clicking on OK button.

III-3 Settings Com




This menu enables the choice of COM number of the serial interface on the PC for the communication with the synthesizer. The COM number can be checked in the system settings of the PC.



The COM number can be changed by clicking on the text field and confirmed by clicking on the OK button.

IV. Subroutine programming



The subroutine is the most important tool for detailed controlling of the synthesis. Every instrument's activity will be defined in subroutines.

The single subroutine consists of 25 steps (columns in the subroutine's table) and every step is defined through the lines: **Time**, **Source**, **Mixed**, **Destin**, **S. Col. Ptr.**, **Delay**, **Branch**. The subroutine's columns define, step by step, from which bottle to which point of the synthesizer for how long the reagent must be delivered.

See the flow diagram of the H-8 synthesizer below:

H-8 - Flow chart

The screenshot shows a software window titled "H-8 on COM_6" with a menu bar (File, Edit, Settings, Control, Window). Inside the window is a table titled "Subroutine Untitled:1.sub". The table has 16 columns and the following rows:

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Time (0.1s)	10	25	5	5	2	6	7	8	9	10						
Source	AMD	TET	ACN	TCA	OXI	CP_A	CP_B	MX_C	GAS	S_1						
Mixed																
Destin	COL	CTOP	MX_B	MX_T	M_W	TRM	TM+T	MX+T	AMD	CL+T						
S. Col. Ptr					ON		ON									
Delay (1 s)			5													
Branch				1					2							

- **Time:** the number in that field defines the duration time of the operation like flow time or wait time. The unit is 0.1 sec. (10 = 1 sec).

- **Source:** the code in the Source field defines from which point the reagent has to be delivered.

Source code table:

Pressed key	Text in Source field	Description
'1'	AMD	Flow from amidite bottle
'2'	TET	Flow from TET bottle
'3'	ACN	Flow from ACN bottle
'4'	TCA	Flow from TCA bottle
'5'	OXI	Flow from OXI bottle
'6'	CP_A	Flow from CP_A bottle
'7'	CP_B	Flow from CP_B bottle
'8'	MX_C	Flow from Mix chamber
'9'	GAS	Flow from ARG connection
'S'	S_1	Flow from U/S_1 bottle
'M'	WMO	Wait only if modifikation
'H'	WTH	Wait only if S-oligo

- **Mixed:** the code in the Mixed field defines from which point the reagent has to be delivered together with the reagent defined in source field. The delivery will be performed by switching the valves for the reagents defined in **both** fields: Source and Mixed. The keys to be pressed and the codes are exactly the same like Source codes.
- **Destin:** the code in the Destin field defines the point to which the reagent has to be delivered.

Source code table:

Pressed key	Text in Destin field	Description
'1'	COL	Flow through the column from the bottom to the waste
'2'	CTOP	Flow through the column from the top to the bottom Works only when GAS as source
'3'	MX_B	Flow through the mix chamber from the bottom to the waste
'4'	MX_T	Flow through the mix chamber from the top to the bottom. Works only when GAS as source
'5'	M_W	Flow through the valve block to the waste.

'6'	TRM	Flow through the column from the bottom to the trityl monitor.
'7'	TM+T	Flow through the column from the bottom to the trityl monitor with flow time correction when more than two columns to deliver simultaneously.
'8'	MX+T	Flow through the mix chamber to the waste. Time correction is not Applicable for this destination.
'9'	CL+T	Flow through the column from the bottom to the waste with flow time correction when more than two columns to deliver simultaneously.
'A'	AMD	Flow into the AMD bottle. Works only when GAS as source

- **S. Col. Ptr.:** means 'Single Column Pointer' which is the pointer for serial execution of the part of the subroutine, which lays between two 'ON' markers in the subroutine. The parts of the subroutine which are beyond the ON markers will be executed simultaneously on all columns. For details see **IV-1-1. Subroutine programming-examples**
- **Delay:** means simply 'wait – do nothing' function. The unit is 1 sec.
- **Branch:** this field contains the branch subroutine's call which will be executed due to the amidite in the synthesis. For details see **IV-1-4. Subroutine programming-examples**
Programming branches

IV-1. Subroutine programming-examples

IV-1-1. Single Column Pointer programming. 1-Col. Mode Parameter.

Assuming all 8 columns synthesis on the H-8 synthesizer and following subroutine:

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
Time (0.1s)	10	10	13	1	7	30									
Source	ACN	GAS	TCA		ACN	GAS									
Mixed															
Destin	M_W	COL	TRM	TRM	COL	COL									
S. Col. Ptr			ON	ON	ON	ON									
Delay (1 s)					3			2							
Branch															

When Parameter **1-Col. Mode** set to **NO**

1-Col. Mode	NO
Switch speed AMD	3
Switch speed Reag.	5
Collect DMT	YES
Add time amd flow	2
Add time reag flow	4
Add time tet/amd flow	2
Add time mix reag flow	3
AMD Dissolve time	0

The subroutine will be executed as follows:

- 1) Step 1: 1 sec ACN flow through the valve block to the waste
- 2) Step 2: 1 sec GAS flow through all 8 columns to the waste (column emptying)

Start of the Single Column Pointer.

- 3) Step 3: 1.3 sec TCA flow through columns 1 & 2 to the trityl monitor.

- 4) Step 4: 0.1 sec wait idle (only for formal marking at the end of the single column pointer).
- 5) Step 3a: 1.3 sec TCA flow through columns 3 & 4 to the trityl monitor.
- 6) Step 4a: 0.1 sec wait idle
- 7) Step 3a: 1.3 sec TCA flow through columns 5 & 6 to the trityl monitor.
- 8) Step 4a: 0.1 sec wait idle
- 9) Step 3a: 1.3 sec TCA flow through columns 7 & 8 to the trityl monitor.
- 10) Step 4a: 0.1 sec wait idle

End of the Single Column Pointer.

- 11) Step 5: 5 sec wait idle.

Start of the Single Column Pointer.

- 12) Step 6: 0.7 sec ACN flow through columns 1 & 2 to the waste.
- 13) Step 7: formal marking at the end of the single column pointer.
- 14) Step 6: 0.7 sec ACN flow through column 3 & 4 to the waste.
- 15) Step 7: formal marking at the end of the single column pointer.
- 16) Step 6: 0.7 sec ACN flow through column 5 & 6 to the waste.
- 17) Step 7: formal marking at the end of the single column pointer.
- 18) Step 6: 0.7 sec ACN flow through column 7 & 8 to the waste.
- 19) Step 7: formal marking at the end of the single column pointer.

End of the Single Column Pointer.

- 20) Step 8: 3 sec GAS flow through all 8 columns to the waste (column emptying)
- 21) Step 8a: 2 sec wait idle.

When Parameter **1-Col. Mode** set to **YES**

1-Col. Mode	YES
Switch speed AMD	3
Switch speed Reag.	5
Collect DMT	YES
Add time amd flow	2
Add time reag flow	4
Add time tet/amd flow	2
Add time mix reag flow	3
AMD Dissolve time	0

- 22) Step 1: 1 sec ACN flow through the valve block to the waste

23) Step 2: 1 sec GAS flow through all 8 columns to the waste (column emptying)

Start of the Single Column Pointer.

- 24) Step 3: 1.3 sec TCA flow through column 1 to the trityl monitor.
- 25) Step 4: 0.1 sec wait idle (only for formal marking of the end of the single column pointer).
- 26) Step 3a: 1.3 sec TCA flow through column 2 the trityl monitor.
- 27) Step 4a: 0.1 sec wait idle
- 28) Step 3a: 1.3 sec TCA flow through column 3 the trityl monitor.
- 29) Step 4a: 0.1 sec wait idle
- 30) Step 3a: 1.3 sec TCA flow through column 4 the trityl monitor.
- 31) Step 4a: 0.1 sec wait idle
- 32) Step 3a: 1.3 sec TCA flow through column 5 the trityl monitor.
- 33) Step 4a: 0.1 sec wait idle
- 34) Step 3a: 1.3 sec TCA flow through column 6 the trityl monitor.
- 35) Step 4a: 0.1 sec wait idle
- 36) Step 3a: 1.3 sec TCA flow through column 7 the trityl monitor.
- 37) Step 4a: 0.1 sec wait idle
- 38) Step 3a: 1.3 sec TCA flow through column 8 the trityl monitor.
- 39) Step 4a: 0.1 sec wait idle

End of the Single Column Pointer.

40) Step 5: 5 sec wait idle.

Start of the Single Column Pointer.

- 41) Step 6: 0.7 sec ACN flow through column 1 to the waste.
- 42) Step 7: formal marking at the end of the single column pointer.
- 43) Step 6: 0.7 sec ACN flow through column 2 to the waste.
- 44) Step 7: formal marking at the end of the single column pointer.
- 45) Step 6: 0.7 sec ACN flow through column 3 to the waste.
- 46) Step 7: formal marking at the end of the single column pointer.
- 47) Step 6: 0.7 sec ACN flow through column 4 to the waste.
- 48) Step 7: formal marking at the end of the single column pointer.
- 49) Step 6: 0.7 sec ACN flow through column 5 to the waste.
- 50) Step 7: formal marking at the end of the single column pointer.
- 51) Step 6: 0.7 sec ACN flow through column 6 to the waste.
- 52) Step 7: formal marking at the end of the single column pointer.
- 53) Step 6: 0.7 sec ACN flow through column 7 to the waste.
- 54) Step 7: formal marking at the end of the single column pointer.
- 55) Step 6: 0.7 sec ACN flow through column 8 to the waste.
- 56) Step 7: formal marking at the end of the single column pointer..

End of the Single Column Pointer.

- 57) Step 8: 3 sec GAS flow through all 8 columns to the waste (column emptying)
- 58) Step 8a: 2 sec wait idle.

The Parameter **1-Col. Mode** set to **YES** means flow only through one column at the same time. This can lead to an increase of cycle times.

IV-1-2. Switch speed AMD-parameter Switch speed Reag.-parameter

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	1
Time (0.1s)	20	10	30	10	20										
Source	CP_A	GAS	AMD	GAS	TET										
Mixed	CP_B		TET		AMD										
Destin	COL	COL	COL	COL	COL										
S. Col. Ptr															
Delay (1 s)	2		2		5										
Branch															

1-Col. Mode	NO
Switch speed AMD	3
Switch speed Reag.	5
Collect DMT	YES
Add time amd flow	2
Add time reag flow	4
Add time tet/amd flow	2
Add time mix reag flow	3
AMD Dissolve time	0

Used when
'Source' or
'Mixed' is AMD

Used when 'Source'
and 'Mixed' not AMD
but any other reagent

The subroutine will be executed as follows:

- 1) Step 1a: flow through all columns to the waste in switching mode: valve CP_A **on**, valve CP_B **off** for 5 units (Switch speed Reag. Parameter)
Step 1b: flow through all columns to the waste in switching mode: valve CP_A **off**, valve CP_B **on** for 5 units.
Repeat Step 1a and 1b until 2 seconds from the field 'Time' are over.
Step 1c: Wait 2 sec. idle.
- 2) Step 2: GAS flow through all columns (emptying of the columns).

- 3) Step 3a: flow through all columns to the waste in switching mode: valve AMD **on**, valve TET **off** for 3 units (Switch speed AMD Parameter)
 Step 3b: flow through all columns to the waste in switching mode: valve AMD **off**, valve TET **on** for 3 units.
 Repeat Step 3a and 3b until 3 seconds from the field 'Time' are over.
 Step 3c: Wait 2 sec. idle.
- 4) Step 4: GAS flow through all columns (emptying of the columns).
- 5) Step 5a: flow through all columns to the waste in switching mode: valve TET **on**, valve AMD **off** for 3 units (Switch speed AMD Parameter)
 Step 5b: flow through all columns to the waste in switching mode: valve TET **off**, valve AMD **on** for 3 units.
 Repeat Step 5a and 5b until 2 seconds from the field 'Time' are over.
 Step 5c: Wait 5 sec. idle.

IV-1-3. Add time amd flow.

Add time reag flow


Add time tet/amd flow

Add time mix reag flow


1-Col. Mode	NO
Switch speed AMD	3
Switch speed Reag.	5
Collect DMT	YES
Add time amd flow	2
Add time reag flow	4
Add time tet/amd flow	2
Add time mix reag flow	3
AMD Dissolve time	0

Assuming all 8 columns synthesis on the H-8 synthesizer and following subroutine:

3) **Step 3:** flow through all columns to the waste for 1 sec + time correction 'add time tet/amd reag flow'. The time correction will be calculated as follows:



4) **Step 4:** flow through all columns to the waste for 1 sec + time correction 'add time mix reag flow'. The time correction will be calculated as follows:



Step 5 & 6: flow to the waste in 'single column mode'. No time correction will be calculated, because maximum 2 columns will be delivered in this mode. See **IV-1-1. Single Column Pointer programming.**

IV-1-4. Branch Programming

'Branch' allows the execution of the subroutine due to the amidite in the sequence and is very useful for coupling programming.

Assuming following sequences in step 6 of the synthesis and subroutine

5'

3'

COL 1:	ACC	AGT	GTA	CCA	CGT	CCG	ACT
COL 2:	CCA	GTG	TAC	CAC	GTA	CGA	CTG
COL 3:	CAG	TGT	ACC	ACG	TAC	GAC	TGG
COL 4:	AGT	GTA	CCA	CGT	ACG	ACT	GGA
COL 5:	GTG	TAC	CAC	GTA	CGA	CTG	GAC
COL 6:	TGT	ACC	ACG	TAC	GAC	TGG	ACA
COL 7:	GTA	CCA	CGT	ACG	ACT	GGA	CAA
COL 8:	TAC	CAC	GTA	CGA	CTG	ACA	AGG

H-8 on COM_6

File Edit Settings Control Window

C:\H-8-MANU\CPL.SUB

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
Time (0.1s)	1	5	10	4					2		10	10	10	15	
Source	TET	ACN	GAS	TET					GAS		ACN	GAS	GAS	ACN	
Mixed															
Destin	M_W	M_W	M_W	COL					CTOP		M_W	M_W	COL	CL+T	
S. Col. Ptr				ON	ON				ON	ON				ON	
Delay (1 s)						2		40							
Branch							2								

When the subroutine execution reaches step 7 of above subroutine:

- 1) The program checks which amidite occurs in the sequence on the column 1. It is **C** in the example.
- 2) The programs check in the Table (III-2 Tables) the subroutine's name from the table on the place pointed by the number from the field 'Branch' in the subroutine.

Method edit

CYCLE			BRANCH TABLE	
Step	Subroutine	Qty	Branch	Subrouti
1	1udet		1	brnch
2	1ucpl		2	1ubrnc
3	1ucap		3	
4	1uoxi		4	
5			5	
6				
7				
8				
9				
10				
11				
12				
13				
14				
15				

Base: C

The number for the branch
From the field 'Branch' in the subroutine

Subroutine's name which will be
executed now

- 3) The program executes completely the subroutine '1ubrnc' which is pointed by the number from the field 'Branch' in the subroutine. The branch subroutine, in this case 1ubrnc, will be executed for **all** sequences in which **C** occurs, in the example for the columns **1 2 5**.

- 4) When branch subroutine is finished, the execution returns to step 7 of the subroutine and checks the next amidite which was not involved yet; **G** on column 3 in the example.
- 5) The programs check in the Table the subroutine's name from the table on the place number 2 for G, which is 1ubrnc.

Method edit

CYCLE			BRANCH TABLE	
Step	Subroutine	Qty		
1	1udet		1ubrnc.sub	Base : G
2	1ucpl		1ucap.sub	Branch Subrouti
3	1ucap		1ucpl.sub	1 brnc
4	1uoxi		1udet.sub	2 1ubrnc
5			1umix.sub	3
6			1uoxi.sub	4
7			brnc.sub	5
8			brnc-i.sub	
9			brncrna.sub	
10			cap.sub	
11			cap-long.sub	
12			caprna.sub	
13			cpl.sub	
14			cpl-long.sub	
15			cplrna.sub	
			det.sub	
			det-long.sub	
			detrna.sub	
			mixed.sub	
			oxi.sub	

- 6) The program executes completely the subroutine '1ubrnc' all sequences in which **G** occurs, in the example for the columns 3 7.
- 7) When branch subroutine is finished, the execution returns to step 7 of the subroutine and checks the next amidite which was not involved yet; **A** on column 4 in the example.
- 8) The programs check in the Table the subroutine's name from the table on the place number 2 for A, which is also 1umix.

Method edit

ROUTINE			BRANCH TABLE	
Step	Subroutine	Qty		
1	1udet		1ubrnc.sub	Base : A
2	1ucpl		1ucap.sub	Branch Subrouti
3	1ucap		1ucpl.sub	1 brnc
4	1uoxi		1udet.sub	2 1umix
5			1umix.sub	3
6			1uoxi.sub	4
7			brnc.sub	5
8			brnc-i.sub	
9			brncrna.sub	
10			cap.sub	
11			cap-long.sub	
12			caprna.sub	
13			cpl.sub	
14			cpl-long.sub	
15			cplrna.sub	
			det.sub	
			det-long.sub	
			detrna.sub	
			mixed.sub	
			oxi.sub	

- 9) The program executes completely the subroutine '1umix' **all** sequences in which **A** occurs, in the example for the columns **4 8**.
- 10) When branch subroutine is finished, the execution returns to step 7 of the subroutine and checks the next amidite which was not involved yet; **T** on column 6 in the example.
- 11) The programs check in the Table the subroutine's name from the table on the place number 2 for T, which is also 1ubrnc.

The screenshot shows the 'Method edit' interface. On the left is a 'CYCLE' table with columns 'Step', 'Subroutine', and 'Qty'. The 'Subroutine' column contains entries for steps 1 through 4: '1udet', '1ucpl', '1ucap', and '1uoxi'. To the right is a list of subroutines, with '1umix.sub' highlighted. On the far right is a 'BRANCH TABLE' with a 'Base' dropdown set to 'T' and a table with columns 'Branch' and 'Subrouti'. The 'BRANCH TABLE' contains entries for branches 1 through 5, with branch 2 pointing to '1ubrnc'. An arrow points from the '1ubrnc' entry in the 'BRANCH TABLE' to the '1umix.sub' entry in the subroutine list.

Step	Subroutine	Qty
1	1udet	
2	1ucpl	
3	1ucap	
4	1uoxi	
5		
6		
7		
8		
9		
10		
11		
12		
13		
14		
15		

Branch	Subrouti
1	brnc
2	1ubrnc
3	
4	
5	

- 12) When branch subroutine is finished, the execution returns to step 7 of the subroutine and checks the next amidite which was not involved yet. The execution resumes at the step 8 of the subroutines when all amidites had been already involved.

It is not allowed to use branch call within branch subroutine.